

Product Datasheet

Recombinant Mouse Rank/Tnfrsf11a/Cd265 Protein, His Tag (orb1551753)

Description	Recombinant Mouse Rank/Tnfrsf11a/Cd265 Protein, His Tag
Species/Host	Mammal
Reactivity	Mouse
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	Tnfrsf11a
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the mouse Tnfrsf11a(Val31-Ser214) was fused with the C-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	O35305
Expiration Date	6 months from date of receipt.