

Product Datasheet

Recombinant Human PTP1C/PTPN6 Protein, His Tag (orb1551765)

Description	Recombinant Human PTP1C/PTPN6 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	PTPN6
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human PTPN6(Lys243-Ile541) was fused with the C-terminal His Tag and expressed in E. coli.
Uniprot ID	P29350
Expiration Date	6 months from date of receipt.