

Product Datasheet

Recombinant Human PTHLH Protein, His Tag (orb1551766)

Description	Recombinant Human PTHLH Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	PTHLH
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human PTHLH(Ala37-Arg175) was fused with the C-terminal His Tag and expressed in E. coli.
Uniprot ID	P12272
Expiration Date	6 months from date of receipt.