

Product Datasheet

Recombinant Mouse Prolactin R/Prlr Protein, Fc Tag (orb1551789)

Description	Recombinant Mouse Prolactin R/Prlr Protein, Fc Tag
Species/Host	Mammal
Reactivity	Mouse
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	Prlr
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the mouse Prlr(Gln20-Asp229) was fused with the C-terminal Fc Tag and expressed in Mammalian cells.
Uniprot ID	Q08501
Expiration Date	6 months from date of receipt.