

Product Datasheet

Recombinant Human PPP1R1A Protein, His Tag (orb1551805)

Description	Recombinant Human PPP1R1A Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	PPP1R1A
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human PPP1R1A(Met1-Gly147) was fused with the C-terminal His Tag and expressed in E. coli.
Uniprot ID	Q13522
Expiration Date	6 months from date of receipt.