

## Product Datasheet

# Recombinant Human PMP2/FABP8 Protein, His Tag (orb1551819)

<b>Description</b>	Recombinant Human PMP2/FABP8 Protein, His Tag
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Human
<b>Endotoxins</b>	Please contact our customer service for more information.
<b>Target</b>	PMP2
<b>Note</b>	For research use only
<b>Purity</b>	≥90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the human PMP2(Met1-Val132) was fused with the N-terminal His Tag and expressed in E. coli.
<b>Uniprot ID</b>	<b>P02689</b>
<b>Expiration Date</b>	6 months from date of receipt.