

Product Datasheet

Recombinant Human PLA2G7/PAF-AH/Lp-PLA2 Protein, His Tag (orb1551827)

Description	Recombinant Human PLA2G7/PAF-AH/Lp-PLA2 Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	PLA2G7
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human PLA2G7(Phe22-Asn441) was fused with the C-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	Q13093
Expiration Date	6 months from date of receipt.