

## Product Datasheet

# Recombinant Human PLA2G16/PLAAT3 Protein, His Tag (orb1551829)

<b>Description</b>	Recombinant Human PLA2G16/PLAAT3 Protein, His Tag
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Human
<b>Endotoxins</b>	Please contact our customer service for more information.
<b>Target</b>	PLAAT3
<b>Note</b>	For research use only
<b>Purity</b>	≥90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the human PLA2G16(Asp12-Asp132) was fused with the N-terminal His Tag and expressed in E. coli.
<b>Uniprot ID</b>	<b>P53816</b>
<b>Expiration Date</b>	6 months from date of receipt.