

Product Datasheet

Recombinant Human PKM2/PKM Protein, His Tag (orb1551830)

Description	Recombinant Human PKM2/PKM Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	PKM
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human PKM(Ser2-Pro531) was fused with the N-terminal His Tag and expressed in E. coli.
Uniprot ID	P14618
Expiration Date	6 months from date of receipt.