

Product Datasheet

Recombinant Human PKI- β Protein, His Tag (orb1551832)

Description	Recombinant Human PKI- β Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	PKIB
Note	For research use only
Purity	$\geq 90\%$ as determined by SDS-PAGE
Source	This protein contains the human PKIB(Met1-Lys78) was fused with the N-terminal His Tag and expressed in E. coli.
Uniprot ID	Q9C010
Expiration Date	6 months from date of receipt.