

## Product Datasheet

# Recombinant *Streptomyces hygroscopicus* Phosphinothricin N-acetyltransferase/bar Protein, No Tag (orb1551843)

<b>Description</b>	Recombinant <i>Streptomyces hygroscopicus</i> Phosphinothricin N-acetyltransferase/bar Protein, No Tag
<b>Species/Host</b>	<i>E. coli</i>
<b>Reactivity</b>	Bacteria
<b>Endotoxins</b>	Please contact our customer service for more information.
<b>Target</b>	PPT N-acetyltransferase
<b>Note</b>	For research use only
<b>Purity</b>	≥90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the <i>Streptomyces hygroscopicus</i> bar(Met1-Ile183) was fused without Tag and expressed in <i>E. coli</i> .
<b>Uniprot ID</b>	<a href="#">P16426</a>
<b>Expiration Date</b>	6 months from date of receipt.