

Product Datasheet

Recombinant Human PFN2 Protein, No Tag (orb1551851)

Description	Recombinant Human PFN2 Protein, No Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	PFN2
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human PFN2(Met1-Phe140) was fused without Tag and expressed in E. coli.
Uniprot ID	P35080
Expiration Date	6 months from date of receipt.