

Product Datasheet

Recombinant Human Peroxiredoxin 5/PRDX5 Protein, His Tag (orb1551857)

Description	Recombinant Human Peroxiredoxin 5/PRDX5 Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	PRDX5
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human PRDX5(Met53-Leu214) was fused with the N-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	P30044
Expiration Date	6 months from date of receipt.