

Product Datasheet

Recombinant Human PRDX4/Peroxiredoxin 4N C-his Protein, His Tag & His Tag (orb1551858)

Description	Recombinant Human PRDX4/Peroxiredoxin 4N C-his Protein, His Tag & His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	PRDX4
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human PRDX4(Met1-Asn271) was fused with the C-terminal His Tag N-terminal His Tag and expressed in E. coli.
Uniprot ID	Q13162
Expiration Date	6 months from date of receipt.