

Product Datasheet

Recombinant Human PDGFR- β /PDGFRB Protein, His Tag (orb1551874)

Description	Recombinant Human PDGFR- β /PDGFRB Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	PDGFRB
Note	For research use only
Purity	$\geq 90\%$ as determined by SDS-PAGE
Source	This protein contains the human PDGFRB(Leu33-Phe530) was fused with the C-terminal His Tag and expressed in E. coli.
Uniprot ID	P09619
Expiration Date	6 months from date of receipt.