

Product Datasheet

Recombinant Human P4HB Protein, His Tag (orb1551909)

Description	Recombinant Human P4HB Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	P4HB
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human P4HB(Asp18-Lys505) was fused with the C-terminal His Tag and expressed in E. coli.
Uniprot ID	P07237
Expiration Date	6 months from date of receipt.