

Product Datasheet

Recombinant Human OX40/TNFRSF4/CD134 Protein, Fc Tag (orb1551913)

Description	Recombinant Human OX40/TNFRSF4/CD134 Protein, Fc Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	TNFRSF4
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human Tnfsf4(Leu29-Ala216) was fused with the N-terminal Fc Tag and expressed in Mammalian cells.
Uniprot ID	P43489
Expiration Date	6 months from date of receipt.