

## Product Datasheet

# Recombinant Mouse OX40 Ligand/Tnfsf4 Protein, His Tag (orb1551918)

<b>Description</b>	Recombinant Mouse OX40 Ligand/Tnfsf4 Protein, His Tag
<b>Species/Host</b>	Mammal
<b>Reactivity</b>	Mouse
<b>Endotoxins</b>	< 1.0 EU per µg of the protein as determined by the LAL method.
<b>Target</b>	Tnfsf4
<b>Note</b>	For research use only
<b>Purity</b>	≥90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the mouse Tnfsf4(Ser51-Leu198) was fused with the N-terminal His Tag and expressed in Mammalian cells.
<b>Uniprot ID</b>	<b>P43488</b>
<b>Expiration Date</b>	6 months from date of receipt.