

Product Datasheet

Recombinant Human NOGO/RTN4 Protein, Gst Tag (orb1551950)

Description	Recombinant Human NOGO/RTN4 Protein, Gst Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	RTN4
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human RTN4(Met1-Val185) was fused with the N-terminal Gst Tag and expressed in E. coli.
Uniprot ID	Q9NQC3
Expiration Date	6 months from date of receipt.