

Product Datasheet

Recombinant Human Neurtin/NRN1 Protein, His Tag (orb1551971)

Description	Recombinant Human Neurtin/NRN1 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	NRN1
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human NRN1(Ala28-Gly116) was fused with the N-terminal His Tag and expressed in E. coli.
Uniprot ID	Q9NPD7
Expiration Date	6 months from date of receipt.