

Product Datasheet

Recombinant Human NECTIN-2/PVRL2/NECTIN1 Protein, His Tag (orb1551977)

Description	Recombinant Human NECTIN-2/PVRL2/NECTIN1 Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	NECTIN2
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human NECTIN2(Gln32-Leu360) was fused with the C-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	Q92692
Expiration Date	6 months from date of receipt.