

Product Datasheet

Recombinant E. coli Mug Protein, No Tag (orb1551996)

Description	Recombinant E. coli Mug Protein, No Tag
Species/Host	E. coli
Reactivity	E. coli
Endotoxins	Please contact our customer service for more information.
Target	mug
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the E. coli. mug(Met1-Arg168) was fused without Tag and expressed in E. coli.
Uniprot ID	P0A9H1
Expiration Date	6 months from date of receipt.