

Product Datasheet

Recombinant Mouse Mrp8/S100a8 Protein, His Tag (orb1552000)

Description	Recombinant Mouse Mrp8/S100a8 Protein, His Tag
Species/Host	E. coli
Reactivity	Mouse
Endotoxins	Please contact our customer service for more information.
Target	S100a8
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the mouse S100a8(Met1-Glu93) was fused with the C-terminal His Tag and expressed in E. coli.
Uniprot ID	P05109
Expiration Date	6 months from date of receipt.