

Product Datasheet

Recombinant Human METAP1D Protein, N/ His Tag (orb1552025)

| | |
|------------------------|---|
| Description | Recombinant Human METAP1D Protein, N/ His Tag |
| Species/Host | E. coli |
| Reactivity | Human |
| Endotoxins | Please contact our customer service for more information. |
| Target | METAP1D |
| Note | For research use only |
| Purity | ≥90% as determined by SDS-PAGE |
| Source | This protein contains the human METAP1D(Arg44-Ala335) was fused with the N/C-terminal His Tag and expressed in E. coli. |
| Uniprot ID | Q6UB28 |
| Expiration Date | 6 months from date of receipt. |