

Product Datasheet

Recombinant Human M-CSF R/CSF1R/CD115 Protein, Fc Tag (orb1552039)

Description	Recombinant Human M-CSF R/CSF1R/CD115 Protein, Fc Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	CSF1R
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human CSF1R(Ile20-Pro517) was fused with the C-terminal Fc Tag and expressed in Mammalian cells.
Uniprot ID	P07333
Expiration Date	6 months from date of receipt.