

Product Datasheet

Recombinant Human MAP1LC3A Protein, His Tag (orb1552048)

Description	Recombinant Human MAP1LC3A Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	MAP1LC3A
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human MAP1LC3A(Met1-Phe121) was fused with the C-terminal His Tag and expressed in E. coli.
Uniprot ID	Q9H492
Expiration Date	6 months from date of receipt.