

Product Datasheet

Recombinant E. coli Lon Protease Protein, His Tag (orb1552077)

Description	Recombinant E. coli Lon Protease Protein, His Tag
Species/Host	E. coli
Reactivity	E. coli
Endotoxins	Please contact our customer service for more information.
Target	lon
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the E. coli. lon(Asn2-Lys784) was fused with the C-terminal His Tag and expressed in E. coli.
Uniprot ID	P0A9M0
Expiration Date	6 months from date of receipt.