

Product Datasheet

Recombinant Human LMCD1 Protein, His Tag (orb1552079)

Description	Recombinant Human LMCD1 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	LMCD1
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human LMCD1(Met1-Ser365) was fused with the C-terminal His Tag and expressed in E. coli.
Uniprot ID	Q9NZU5
Expiration Date	6 months from date of receipt.