

## Product Datasheet

# Recombinant Human LIF/Leukemia inhibitory factor Protein, No Tag (orb1552092)

<b>Description</b>	Recombinant Human LIF/Leukemia inhibitory factor Protein, No Tag
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Human
<b>Endotoxins</b>	Please contact our customer service for more information.
<b>Target</b>	LIF
<b>Note</b>	For research use only
<b>Purity</b>	≥90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the human LIF(Ser23-Phe202) was fused without Tag and expressed in E. coli.
<b>Uniprot ID</b>	<b>P15018</b>
<b>Expiration Date</b>	6 months from date of receipt.