

Product Datasheet

Recombinant Mouse Lif Protein, No Tag (orb1552093)

Description	Recombinant Mouse Lif Protein, No Tag
Species/Host	E. coli
Reactivity	Mouse
Endotoxins	Please contact our customer service for more information.
Target	Lif
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the mouse Lif(Ser24-Phe203) was fused without Tag and expressed in E. coli.
Uniprot ID	P09056
Expiration Date	6 months from date of receipt.