

Product Datasheet

Recombinant Human Leptin/LEP Protein, No Tag (orb1552095)

Description	Recombinant Human Leptin/LEP Protein, No Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	LEP
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human LEP(Val22-Cys167) was fused without Tag and expressed in E. coli.
Uniprot ID	P41159
Expiration Date	6 months from date of receipt.