

Product Datasheet

Recombinant Human LCP2 Protein, His Tag (orb1552105)

| | |
|------------------------|---|
| Description | Recombinant Human LCP2 Protein, His Tag |
| Species/Host | E. coli |
| Reactivity | Human |
| Endotoxins | Please contact our customer service for more information. |
| Target | LCP2 |
| Note | For research use only |
| Purity | ≥90% as determined by SDS-PAGE |
| Source | This protein contains the human LCP2(Met1-Pro533) was fused with the C-terminal His Tag and expressed in E. coli. |
| Uniprot ID | Q13094 |
| Expiration Date | 6 months from date of receipt. |