

Product Datasheet

Recombinant E. coli LacI Protein, His Tag (orb1552122)

Description	Recombinant E. coli LacI Protein, His Tag
Species/Host	E. coli
Reactivity	E. coli
Endotoxins	Please contact our customer service for more information.
Target	lacI
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the E. coli. lacI(Met1-Gln360) was fused with the C-terminal His Tag and expressed in E. coli.
Uniprot ID	J7QMK7
Expiration Date	6 months from date of receipt.