

Product Datasheet

Recombinant Human Kallikrein 11/KLK11 Protein, His Tag (orb1552144)

Description	Recombinant Human Kallikrein 11/KLK11 Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	KLK11
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human KLK11(Glu19-Asn250) was fused with the C-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	Q9UBX7
Expiration Date	6 months from date of receipt.