

Product Datasheet

Recombinant Human ITPase/ITPA Protein, His Tag (orb1552152)

| | |
|------------------------|-------------------------------------------------------------------------------------------------------------------|
| Description | Recombinant Human ITPase/ITPA Protein, His Tag |
| Species/Host | E. coli |
| Reactivity | Human |
| Endotoxins | Please contact our customer service for more information. |
| Target | ITPA |
| Note | For research use only |
| Purity | ≥90% as determined by SDS-PAGE |
| Source | This protein contains the human ITPA(Ala2-Ala194) was fused with the C-terminal His Tag and expressed in E. coli. |
| Uniprot ID | Q9BY32 |
| Expiration Date | 6 months from date of receipt. |