

## Product Datasheet

# Recombinant Human INHBC Protein, His Tag (orb1552167)

<b>Description</b>	Recombinant Human INHBC Protein, His Tag
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Human
<b>Endotoxins</b>	Please contact our customer service for more information.
<b>Target</b>	INHBC
<b>Note</b>	For research use only
<b>Purity</b>	≥90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the human INHBC(Thr19-Ser352) was fused with the C-terminal His Tag and expressed in E. coli.
<b>Uniprot ID</b>	<b>P55103</b>
<b>Expiration Date</b>	6 months from date of receipt.