

Product Datasheet

Recombinant Human IL-31 Protein, No Tag (orb1552175)

Description	Recombinant Human IL-31 Protein, No Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	IL31
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human IL31(Ser24-Thr164) was fused without Tag and expressed in E. coli.
Uniprot ID	Q6EBC2
Expiration Date	6 months from date of receipt.