

## Product Datasheet

# Recombinant Mouse Il-23/Il23a Protein, No Tag (orb1552176)

<b>Description</b>	Recombinant Mouse Il-23/Il23a Protein, No Tag
<b>Species/Host</b>	Mammal
<b>Reactivity</b>	Mouse
<b>Endotoxins</b>	< 1.0 EU per µg of the protein as determined by the LAL method.
<b>Target</b>	Il23a
<b>Note</b>	For research use only
<b>Purity</b>	≥90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the mouse Il23a(Val22-Ala196&Met23-Ser335) was fused without Tag and expressed in Mammalian cells.
<b>Uniprot ID</b>	<a href="#">Q9EQ14</a>
<b>Expiration Date</b>	6 months from date of receipt.