

Product Datasheet

Recombinant Mouse IL-1rn Protein, No Tag (orb1552181)

Description	Recombinant Mouse IL-1rn Protein, No Tag
Species/Host	E. coli
Reactivity	Mouse
Endotoxins	Please contact our customer service for more information.
Target	Il1rn
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the mouse Il1rn(Arg27-Gln178) was fused without Tag and expressed in E. coli.
Uniprot ID	P25085
Expiration Date	6 months from date of receipt.