

## Product Datasheet

# Recombinant Mouse Il-12B/Il-12p40/Il12b Protein, His Tag (orb1552191)

<b>Description</b>	Recombinant Mouse Il-12B/Il-12p40/Il12b Protein, His Tag
<b>Species/Host</b>	Mammal
<b>Reactivity</b>	Mouse
<b>Endotoxins</b>	< 1.0 EU per µg of the protein as determined by the LAL method.
<b>Target</b>	Il12b
<b>Note</b>	For research use only
<b>Purity</b>	≥90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the mouse Il12b(Met23-Ser335) was fused with the C-terminal His Tag and expressed in Mammalian cells.
<b>Uniprot ID</b>	<a href="#">P43432</a>
<b>Expiration Date</b>	6 months from date of receipt.