

## Product Datasheet

# Recombinant Mouse Il-12/Il12a Protein, No Tag (orb1552192)

<b>Description</b>	Recombinant Mouse Il-12/Il12a Protein, No Tag
<b>Species/Host</b>	Mammal
<b>Reactivity</b>	Mouse
<b>Endotoxins</b>	< 1.0 EU per µg of the protein as determined by the LAL method.
<b>Target</b>	Il12a
<b>Note</b>	For research use only
<b>Purity</b>	≥90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the mouse Il12a(Met23-Ser335&Arg23-Ala215) was fused without Tag and expressed in Mammalian cells.
<b>Uniprot ID</b>	<a href="#">P43431</a>
<b>Expiration Date</b>	6 months from date of receipt.