

Product Datasheet

Recombinant Human IGF-I/IGF1 Protein, No Tag (orb1552205)

Description	Recombinant Human IGF-I/IGF1 Protein, No Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	IGF1
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human IGF1(Gly49-Ala118) was fused without Tag and expressed in E. coli.
Uniprot ID	P05019
Expiration Date	6 months from date of receipt.