

## Product Datasheet

# Recombinant Human IGF2BP2 Protein, T7 Tag & His Tag (orb1552209)

<b>Description</b>	Recombinant Human IGF2BP2 Protein, T7 Tag & His Tag
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Human
<b>Endotoxins</b>	Please contact our customer service for more information.
<b>Target</b>	IGF2BP2
<b>Note</b>	For research use only
<b>Purity</b>	≥90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the human IGF2BP2(Met1-Thr220) was fused with the N-terminal T7 Tag, C-terminal His Tag and expressed in E. coli.
<b>Uniprot ID</b>	<b>Q9Y6M1</b>
<b>Expiration Date</b>	6 months from date of receipt.