

## Product Datasheet

# Recombinant Human IFN-gamma/IFN- $\gamma$ /IFNG Protein, No Tag (orb1552216)

<b>Description</b>	Recombinant Human IFN-gamma/IFN- $\gamma$ /IFNG Protein, No Tag
<b>Species/Host</b>	Mammal
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1.0 EU per $\mu$ g of the protein as determined by the LAL method.
<b>Target</b>	IFNG
<b>Note</b>	For research use only
<b>Purity</b>	$\geq$ 90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the human IFNG(Gln24-Gln166) was fused without Tag and expressed in Mammalian cells.
<b>Uniprot ID</b>	<b>P01579</b>
<b>Expiration Date</b>	6 months from date of receipt.