

## Product Datasheet

# Recombinant Mouse IFN-gamma/IFN- $\gamma$ /Ifng Protein, No Tag (orb1552217)

<b>Description</b>	Recombinant Mouse IFN-gamma/IFN- $\gamma$ /Ifng Protein, No Tag
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Mouse
<b>Endotoxins</b>	Please contact our customer service for more information.
<b>Target</b>	Ifng
<b>Note</b>	For research use only
<b>Purity</b>	$\geq$ 90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the mouse Ifng(His23-Cys155) was fused without Tag and expressed in E. coli.
<b>Uniprot ID</b>	<b>P01580</b>
<b>Expiration Date</b>	6 months from date of receipt.