

Product Datasheet

Recombinant Human IFN-gamma/IFN- γ /IFNG Protein, No Tag (orb1552218)

Description	Recombinant Human IFN-gamma/IFN- γ /IFNG Protein, No Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	IFNG
Note	For research use only
Purity	$\geq 90\%$ as determined by SDS-PAGE
Source	This protein contains the human IFNG(Gln24-Gln166) was fused without Tag and expressed in E. coli.
Uniprot ID	P01579
Expiration Date	6 months from date of receipt.