

Product Datasheet

Recombinant Mouse IFN-gamma/IFN- γ /Ifng Protein, His Tag (orb1552219)

Description	Recombinant Mouse IFN-gamma/IFN- γ /Ifng Protein, His Tag
Species/Host	Mammal
Reactivity	Mouse
Endotoxins	< 1.0 EU per μ g of the protein as determined by the LAL method.
Target	Ifng
Note	For research use only
Purity	\geq 90% as determined by SDS-PAGE
Source	This protein contains the mouse Ifng(His23-Cys155) was fused with the C-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	P01580
Expiration Date	6 months from date of receipt.