

## Product Datasheet

# Recombinant Human IFN-alpha/beta R2/IFNAR2 Protein, Fc Tag (orb1552224)

<b>Description</b>	Recombinant Human IFN-alpha/beta R2/IFNAR2 Protein, Fc Tag
<b>Species/Host</b>	Mammal
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1.0 EU per µg of the protein as determined by the LAL method.
<b>Target</b>	IFNAR2
<b>Note</b>	For research use only
<b>Purity</b>	≥90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the human IFNAR2(Ile27-Lys243) was fused with the C-terminal Fc Tag and expressed in Mammalian cells.
<b>Uniprot ID</b>	<a href="#">P48551</a>
<b>Expiration Date</b>	6 months from date of receipt.