

Product Datasheet

Recombinant Mouse Ido/Ido1 Protein, His Tag (orb1552229)

Description	Recombinant Mouse Ido/Ido1 Protein, His Tag
Species/Host	E. coli
Reactivity	Mouse
Endotoxins	Please contact our customer service for more information.
Target	Ido1
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the mouse Ido1(Met1-Pro407) was fused with the N-terminal His Tag and expressed in E. coli.
Uniprot ID	P28776
Expiration Date	6 months from date of receipt.