

Product Datasheet

Recombinant Mouse Hvem/Tnfrsf14 Protein, Fc Tag (orb1552243)

Description	Recombinant Mouse Hvem/Tnfrsf14 Protein, Fc Tag
Species/Host	Mammal
Reactivity	Mouse
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	Tnfrsf14
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the mouse Tnfrsf14(Gln39-Val207) was fused with the C-terminal Fc Tag and expressed in Mammalian cells.
Uniprot ID	Q71F55
Expiration Date	6 months from date of receipt.